



Very low viscosity addition curing silicone

For professional use in dental laboratory

Usage: Indicated for duplication in dental laboratory

DUPLIMED

- pink colour
- hardness 14 Sh.A

DUPLIMED HARD

- orange colour
- hardness 25 Sh.A

Characteristics:

- rapid polymerization
- high fluidity
- high index of flowability
- excellent elastic recovery
- excellent reproduction of details
- high dimensional stability
- good resistance to tearing
- compatibility with all technical materials : plaster, coatings and resins

Packaging in flask kg. 1 + 1 and bulk kg. 5 + 5

Mixing ratio 1:1