



Addition cured vinylpolysiloxane for the reproduction of gingival masks in technical laboratory Rigid, be trimmed, scannable type

Usage: Indicated for implant prosthetics

- for direct method
- realization of gingival reproduction in pink vinylpolysiloxane
- easy to trim, for the determination of the shape of the cervical contours and the emergence profile in fixed or mobile prosthesis

Characteristic:

- trimming after 10 minutes from extrusion
- high index of flowability
- opaque colour which give it a natural appearance like real gingiva
- perfectly reflective surface, suitable for scanning with CAD/CAM systems
- balanced tear strenght
- hardness 75 Sh. A after 30 minutes
- high dimensional stability
- no greasy surface

Packaging: cartridge 50 ml.

Extrusion with appropriate dispenser