



## Rigid C- Silicone

For professional use in Dental Laboratory

### Applications:

- Realization of masks in dental laboratory
- Matrix support for the reproduction of gingival tissue (e.g. Meditaly's Gengilla)
- Matrix for use in making temporary bridges/crowns

### Advantages:

- Highly millable without crumbling
- Excellent flowability
- No deformation
- Heat resistance 110°C.
- No greasy surface

### Characteristics:

- Final Hardness 75 Sh. A
- Natural colour
- Mixing with gel or liquid catalyst

**Standard Packaging: 1 jar 5 Kg. + 2 Catalyst gel 60 gr. / 2 Catalyst liquid 30 ml.**