



**Very Low viscosity addition silicone
with high resistance to tearing**
For professional use in dental laboratory

SAR1012

- colour blue
- hardness 12 ShA

SAR 1022

- colour light blue
- hardness 22 ShA

SAR1027

- colour green
- hardness 27 ShA

Usage: Indicated for duplication with high undercuts in dental laboratory

Characteristics:

- rapid polymerization
- high fluidity
- high index of flowability
- excellent elastic recovery
- excellent reproduction of details
- high dimensional stability
- excellent resistance to tearing
- compatibility with all technical materials : plaster, coatings and resins

Packaging in flask kg. 1 + 1 and bulk kg. 5 + 5

Mixing ratio 1:1