



Low viscosity hydrophilic addition silicone (vinylpolysiloxane) for fixed protheses

For professional use in dental clinic

Usage: Most suitable for disclosing areas of poor fit of protheses. Recommended for checking pressure spots and the accuracy of fit of crowns and bridges.

Characteristic:

- milky white opaque
- high index of flowability
- quick, easy accurate auto-mixing application
- minimum waste of material thanks to mini mixing tips
- in case of few elements, possibility of manual mixing
- optimized mixing time
- doesn't create thickness
- non-toxic
- white colour

Packaging: cartridge 50 ml.

Mixing ratio 1:1