



CE

Very high viscosity A-Silicone for fixed protheses

For professional use in dental clinic

Usage: Recommended as special first impression material in case of undercuts

- for the two-stage technique in conjunction with a “ wash” like Extreme Light Hydro or for the single-stage technique (simultaneous) in conjunction with a “ wash” like Extreme Regular Hydro

Characteristic:

- high elasticity and tear strenght
- high index of flowability
- doesn't compress the gingival mucosa
- excellent elastic recovery
- high resistance to deformation
- final hardness 52 Sh. A
- high dimensional stability
- no greasy surface
- violet colour

Packaging: 2 jars (Catalyst + Base) 300+300 ml.