



Medium Viscosity hydrophilic addition silicone

(vinylpolysiloxane) for fixed protheses

For professional use in dental clinic

Usage: Recommended for precision impression as correction fluid

- suitable for the two-stage technique (double impression) or as a monophasic material in conjunction with a putty like Extreme Super Soft Normal

Characteristic:

- fast polymerization for best patient comfort
- high index of flowability
- excellent reproduction of the details
- excellent elastic recovery
- high resistance to deformation
- high dimensional stability
- high tear strength
- hydrophilic - thixotropic
- non-toxic - disinfectable
- bordeaux colour

Packaging : cartridge 50 ml.

Mixing ratio 1:1