



## Very high viscosity A-Silicone for fixed protheses



For professional use in dental clinic

**Usage:** Recommended as first impression material

- in the two-stage technique ( double impression )

### Characteristics:

- rapid polymerization, shorter time in the oral cavity for a greater patient comfort
- high index of flowability
- doesn't compress the gingival mucosa
- excellent elastic recovery
- high resistance to deformation
- final hardness 52 Sh. A
- high dimensional stability
- no greasy surface
- yellow colour

Packaging : 2 jars ( Catalyst + Base ) 300+300 ml.