



Very high viscosity A- Silicone for fixed protheses.

CE

For professional use in dental clinic

Usage: Recommended as first impression material

- for all impression techniques : two-stage (double impression) or single-stage (simultaneous) in conjunction with a “ wash” as Extreme Regular Hydro

Characteristic:

- high index of flowability
- doesn't compress the gingival mucosa
- excellent elastic recovery
- high resistance to deformation
- Final hardness 62 Sh. A
- high dimensional stability
- No greasy surface
- green colour

Packaging : 2 jars (Catalyst + Base) 300+300 ml.