



Transparent A- Silicones

KRISTAL SIL FLUID CE

Indication: Indicated in the realization masks for the moulding of light-curing materials

- in technique with use of transparent flask thanks to a very fluid and flowing material.

Advantages:

- flowability calibrated for use in flask
- high index of fluency
- very high transparency level
- compact and smooth surfaces
- calibrated tear strenght
- high final hardness (70 Sh.A)
- dimensionally stable in time

Packaging: cartridge 50 ml.

Extrusion with appropriate dispenser

KRISTAL SIL TIXO CE

Indication: Indicated in the realization masks for the moulding of light-curing materials

- suitable in free-hand technique in dentistry and dental laboratory thanks to its high thixotropy

Advantages:

- excellent thixotropy
- high index of fluency
- very high transparency level
- compact and smooth surfaces
- calibrated tear strenght
- high final hardness (70 Sh.A)
- dimensionally stable in time